



Production of hair care products



Human hair are very sensitive to the pH variations of applied products. Apart the safety reasons, the pH range also has a big impact on performance of rinse-off cosmetic products. The control of pH is obligatory test in shampoos production. For these purposes, the CHS ChemClean Lab® pH electrode is used. It is operating in wide range of temperatures and pH.

Industry	Hair care cosmetic production
Product	Shampoo



Laboratory pH measurement



Laboratory pH measurement

pH electrode	CHS ChemClean Lab®, S7, part number CHS-5671-NNS
Benefits	<ul style="list-style-type: none">• The color polymer allows easy visual inspection of the electrode• Maintenance-free combined electrode• Suitable for measuring pH at low temperatures• Durable ring PTFE diaphragm minimizes the possibility of clogging and thus interruption of measurement