



The table below lists the specifications for all of SKC's UKEX, ATEX, ANZEx and IECEx certified products for use in the UK, Europe, Australia / New Zealand and Internationally. Please review the product specifications and ensure that your chosen product meets the requirements for use in your Ex potentially explosive atmosphere hazardous location.

Product Range	Models	Certificate(s)	Product marking	ATEX Directive 2014/34/EU		EN IEC 60079-0						EN IEC 60079-10-1 EN 60079-10-2
				Equipment Group	Equipment Category	Level of Protection	Group	Temperature Class	EPL	Ambient Temperature Range, Ta	IP Rating	Zones
Air Sampling Pumps												
Sidekick	224-51MTX 224-52MTX	Sira 02 ATEX 2135	Ⓔ II 1G Ex ia IIC T4 Ga	II	1G	ia (intrinsic safety)	IIA, IIB & IIC	T1, T2, T3 & T4	Ga	-20°C to +40°C	IP 20	0, 1 & 2
Universal	224-44MTX 224-PCMTX4 224-PCMTX8	Sira 02 ATEX 2157	Ⓔ II 1G Ex ia IIC T4 Ga	II	1G	ia (intrinsic safety)	IIA, IIB & IIC	T1, T2, T3 & T4	Ga	-20°C to +40°C	IP 20	0, 1 & 2
AirChek 3000	210-3311	Sira 11 ATEX 2121X IECEX SIR 11.0055X	Ⓔ II 1D Ex ia IIIC T120°C Da	II	1D	ia (intrinsic safety)	IIIA, IIIB & IIIC	T120°C	Da	0°C to +40°C	IP 64	20, 21 & 22
			Ⓔ II 1G Ex ia IIC T4 Ga	II	1G	ia (intrinsic safety)	IIA, IIB & IIC	T1, T2, T3 & T4	Ga	0°C to +40°C	IP 64	0, 1 & 2
			Ⓔ I M1 Ex ia I Ma	I	M1	ia (intrinsic safety)	I	T1, T2, T3 & T4	Ma	0°C to +40°C	IP 64	N/A
	210-3311AZ	ANZEx 12.2007X	Ex ia IIIC T120°C Da	N/A	N/A	ia (intrinsic safety)	IIIA, IIIB & IIIC	T120°C	Da	0°C to +40°C	IP 64	20, 21 & 22
			Ex ia IIC T4 Ga	N/A	N/A	ia (intrinsic safety)	IIA, IIB & IIC	T1, T2, T3 & T4	Ga	0°C to +40°C	IP 64	0, 1 & 2
			Ex ia I Ma	N/A	N/A	ia (intrinsic safety)	I	T1, T2, T3 & T4	Ma	0°C to +40°C	IP 64	N/A
Pocket Pump Touch	220-1000TC	DEMKO 16 ATEX 1600 IECEX UL 16.0113	Ⓔ II 1G Ex ia IIC T4 Ga	II	1G	ia (intrinsic safety)	IIA, IIB & IIC	T1, T2, T3 & T4	Ga	-20°C to +45°C	N/A	0, 1 & 2
AirChek Essential	220-3000	DEMKO 16 ATEX 2288 IECEX UL 19.0100	Ⓔ II 1G Ex ia IIC T4 Ga	II	1G	ia (intrinsic safety)	IIA, IIB & IIC	T1, T2, T3 & T4	Ga	-20°C to +45°C	IP 6X	0, 1 & 2
AirChek Connect	220-4000	DEMKO 16 ATEX 2288 IECEX UL 19.0100	Ⓔ II 1G Ex ia IIC T4 Ga	II	1G	ia (intrinsic safety)	IIA, IIB & IIC	T1, T2, T3 & T4	Ga	-20°C to +45°C	IP 6X	0, 1 & 2



Product Range	Models	Certificate(s)	Product marking	ATEX Directive 2014/34/EU		EN IEC 60079-0						EN IEC 60079-10-1 EN 60079-10-2	
				Equipment Group	Equipment Category	Level of Protection	Group	Temperature Class	EPL	Ambient Temperature Range, Ta	IP Rating	Zones	
Noise Dosimeters													
NoiseChek	701-001S 701-001NBS	UL 21 ATEX 2537 IECEX UL 21.0026	 II 1G Ex ia IIC T4 Ga	II	1G	ia (intrinsic safety)	IIA, IIB & IIC	T1, T2, T3 & T4	Ga	-20°C to +45°C	IP 6X	0, 1 & 2	